

Dengue Virus NS1 Polyclonal Antibody HRP Conjugated



Catalog No: #C91040HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dengue Virus NS1 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC
Species Reactivity	Dengue Virus 2
Immunogen Type	801-900/3391
Conjugates	HRP
Target Name	Dengue Virus NS1
Other Names	Dengue Virus NS1
Accession No.	Uniprot:P29990Gene ID:1494449
Uniprot	P29990
GeneID	1494449
Calculated MW	40
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only