

DCUN1D1 Polyclonal Antibody FITC Conjugated

Catalog No: #C91043FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DCUN1D1 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms
Immunogen Type	1-100/259
Conjugates	FITC
Target Name	DCUN1D1
Other Names	DCUN1D1
Accession No.	Uniprot:Q96GG9 Gene ID:54165
Uniprot	Q96GG9
GeneID	54165
Calculated MW	30
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only