

CD68 Polyclonal Antibody RBITC Conjugated

Catalog No: #C91045RBITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD68 Polyclonal Antibody RBITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu
Immunogen Type	1-130/354
Conjugates	RBITC
Target Name	CD68
Other Names	CD68
Accession No.	Uniprot:P34810Gene ID:968
Uniprot	P34810
GeneID	968
Calculated MW	37
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only