

Cytokeratin 5 Polyclonal Antibody PE Conjugated

Catalog No: #C91049PE



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Cytokeratin 5 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Immunogen Type	451-580/580
Conjugates	PE
Target Name	Cytokeratin 5
Other Names	Cytokeratin 5
Accession No.	Uniprot:P13647Gene ID:3852
Uniprot	P13647
GeneID	3852
Calculated MW	64
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only