

Jagged1 Polyclonal Antibody Cy5 Conjugated

Catalog No: #C91054Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Jagged1 Polyclonal Antibody Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1101-1218/1218
Conjugates	Cy5
Target Name	Jagged1
Other Names	Jagged1
Accession No.	Uniprot:P78504Gene ID:182
Uniprot	P78504
GeneID	182
Calculated MW	134
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only