

EFR3B Polyclonal Antibody HRP Conjugated

Catalog No: #C91056HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EFR3B Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	21-120/817
Conjugates	HRP
Target Name	EFR3B
Other Names	EFR3B
Accession No.	Uniprot:Q9Y2G0Gene ID:22979
Uniprot	Q9Y2G0
GeneID	22979
Calculated MW	92
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only