

EHD3 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91057Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EHD3 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	11-110/535
Conjugates	Bio
Target Name	EHD3
Other Names	EHD3
Accession No.	Uniprot:Q9NZN3Gene ID:30845
Uniprot	Q9NZN3
GeneID	30845
Calculated MW	61
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only