

## IL-23R Polyclonal Antibody PE-Cy5 Conjugated

Catalog No: #C91066PE-Cy5

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	IL-23R Polyclonal Antibody PE-Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	401-550/629
Conjugates	PE-Cy5
Target Name	IL-23R
Other Names	IL-23R
Accession No.	Uniprot:Q5VWK5Gene ID:149233
Uniprot	Q5VWK5
GeneID	149233
Calculated MW	67
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only