

ERRF11 Polyclonal Antibody Cy5 Conjugated

Catalog No: #C91068Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ERRF11 Polyclonal Antibody Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/462
Conjugates	Cy5
Target Name	ERRF11
Other Names	ERRF11
Accession No.	Uniprot:Q9UJM3Gene ID:54206
Uniprot	Q9UJM3
GeneID	54206
Calculated MW	51
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only