

Exportin T Polyclonal Antibody Biotin Conjugated

Catalog No: #C91070Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Exportin T Polyclonal Antibody Biotin Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | WB IHC ICC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 1-100/962 |
| Conjugates | Bio |
| Target Name | Exportin T |
| Other Names | Exportin T |
| Accession No. | Uniprot:O43592Gene ID:11260 |
| Uniprot | O43592 |
| GeneID | 11260 |
| Calculated MW | 110 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only