

MDR1/P Glycoprotein Polyclonal Antibody Cy5 Conjugated

Catalog No: #C91072Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MDR1/P Glycoprotein Polyclonal Antibody Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1051-1280/1280
Conjugates	Cy5
Target Name	MDR1/P Glycoprotein
Other Names	MDR1/P Glycoprotein
Accession No.	Uniprot:P08183Gene ID:5243
Uniprot	P08183
GeneID	5243
Calculated MW	141
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only