

# PAFR/PAF Receptor Polyclonal Antibody APC Conjugated

Catalog No: #C91077APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PAFR/PAF Receptor Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	231-342/342
Conjugates	APC
Target Name	PAFR/PAF Receptor
Other Names	PAFR/PAF Receptor
Accession No.	Uniprot:P25105Gene ID:5724
Uniprot	P25105
GeneID	5724
Calculated MW	38
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only