

MUC1 Polyclonal Antibody PE-Cy5.5 Conjugated

Catalog No: #C91084PE-Cy5.5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MUC1 Polyclonal Antibody PE-Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1101-1255/1255
Conjugates	PE-Cy5.5
Target Name	MUC1
Other Names	MUC1
Accession No.	Uniprot:P15941Gene ID:4582
Uniprot	P15941
GeneID	4582
Calculated MW	20/138
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only