

CXCL10/IP10 Polyclonal Antibody PE-Cy5.5 Conjugated

Catalog No: #C91090PE-Cy5.5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CXCL10/IP10 Polyclonal Antibody PE-Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	35-98/98
Conjugates	PE-Cy5.5
Target Name	CXCL10/IP10
Other Names	CXCL10/IP10
Accession No.	Uniprot:P17515Gene ID:15945
Uniprot	P17515
GeneID	15945
Calculated MW	10
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only