

CD7 Polyclonal Antibody PE Conjugated

Catalog No: #C91105PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD7 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Conjugates	PE
Target Name	CD7
Other Names	CD7
Accession No.	Uniprot:P09564Gene ID:924
Uniprot	P09564
GeneID	924
Calculated MW	24
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only