

phospho-ERK1 / 2 (Thr183/Tyr185) Polyclonal Antibody Biotin Conjugated

Catalog No: #C91108Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-ERK1 / 2 (Thr183/Tyr185) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	FL(p-T)E(p-Y)V
Conjugates	Bio
Target Name	phospho-ERK1 / 2 (Thr183/Tyr185)
Other Names	phospho-ERK1 / 2 (Thr183/Tyr185)
Accession No.	Uniprot:P63085Gene ID:5595
Uniprot	P63085
GeneID	5595
Calculated MW	42/44
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only