

GPCR EX33 Polyclonal Antibody PE Conjugated

Catalog No: #C91115PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPCR EX33 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms
Immunogen Type	21-120/369
Conjugates	PE
Target Name	GPCR EX33
Other Names	GPCR EX33
Accession No.	Uniprot:Q9NQS5Gene ID:80910
Uniprot	Q9NQS5
GeneID	80910
Calculated MW	44
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only