

Alkaline Phosphatase, Tissue Non-Specific isozyme Polyclonal Antibody PE Conjugated



Catalog No: #C91117PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alkaline Phosphatase, Tissue Non-Specific isozyme Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	52-150/524
Conjugates	PE
Target Name	Alkaline Phosphatase, Tissue Non-Specific isozyme
Other Names	Alkaline Phosphatase, Tissue Non-Specific isozyme
Accession No.	Uniprot:P09242Gene ID:249
Uniprot	P09242
GeneID	249
Calculated MW	56
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only