

GPR177 Polyclonal Antibody PE Conjugated

Catalog No: #C91122PE



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	GPR177 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	351-450/541
Conjugates	PE
Target Name	GPR177
Other Names	GPR177
Accession No.	Uniprot:Q5T9L3Gene ID:79971
Uniprot	Q5T9L3
GeneID	79971
Calculated MW	62
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only