

IL-7Ra Polyclonal Antibody PE Conjugated

Catalog No: #C91128PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-7Ra Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	251-350/459
Conjugates	PE
Target Name	IL-7Ra
Other Names	IL-7Ra
Accession No.	Uniprot:P16871Gene ID:3575
Uniprot	P16871
GeneID	3575
Calculated MW	50
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only