

# HBsAg Polyclonal Antibody bs-1557G-Biotin Conjugated

Catalog No: #C91155bs-1557G-Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HBsAg Polyclonal Antibody bs-1557G-Biotin Conjugated
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu
Conjugates	bs-1557G-Bio
Target Name	HBsAg
Other Names	HBsAg
Accession No.	Uniprot:Q67926Gene ID:944569
Uniprot	Q67926
GeneID	944569
Calculated MW	44
Formulation	0.01M PBS(pH7.4) with 0.03% Proclin300.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only