

GLP-1R Polyclonal Antibody RBITC Conjugated

Catalog No: #C91159RBITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GLP-1R Polyclonal Antibody RBITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/463
Conjugates	RBITC
Target Name	GLP-1R
Other Names	GLP-1R
Accession No.	Uniprot:P32301Gene ID:25051
Uniprot	P32301
GeneID	25051
Calculated MW	51
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only