

SFTPD Polyclonal Antibody FITC Conjugated

Catalog No: #C91168FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SFTPD Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	301-375/375
Conjugates	FITC
Target Name	SFTPD
Other Names	SFTPD
Accession No.	Uniprot:P35247Gene ID:6441
Uniprot	P35247
GeneID	6441
Calculated MW	41
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only