

FAM83H Polyclonal Antibody FITC Conjugated

Catalog No: #C91172FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FAM83H Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1001-1179/1179
Conjugates	FITC
Target Name	FAM83H
Other Names	FAM83H
Accession No.	Uniprot:Q6ZRV2Gene ID:286077
Uniprot	Q6ZRV2
GeneID	286077
Calculated MW	127
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only