

ABCB5 Polyclonal Antibody FITC Conjugated

Catalog No: #C91173FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ABCB5 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	841-883/1257
Conjugates	FITC
Target Name	ABCB5
Other Names	ABCB5
Accession No.	Uniprot:Q2M3G0Gene ID:340273
Uniprot	Q2M3G0
GeneID	340273
Calculated MW	89/138
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:200-800

Note: This product is for in vitro research use only