

FGD3 Polyclonal Antibody FITC Conjugated

Catalog No: #C91175FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FGD3 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	501-600/725
Conjugates	FITC
Target Name	FGD3
Other Names	FGD3
Accession No.	Uniprot:Q5JSP0Gene ID:89846
Uniprot	Q5JSP0
GeneID	89846
Calculated MW	79
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only