

Cadherin28 Polyclonal Antibody FITC Conjugated

Catalog No: #C91179FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cadherin28 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	FITC
Target Name	Cadherin28
Other Names	Cadherin28
Accession No.	Uniprot:Q6ZTQ4Gene ID:222256
Uniprot	Q6ZTQ4
GeneID	222256
Calculated MW	96
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only