

RUFY4/FLJ46536 Polyclonal Antibody FITC Conjugated

Catalog No: #C91182FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RUFY4/FLJ46536 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	131-230/372
Conjugates	FITC
Target Name	RUFY4/FLJ46536
Other Names	RUFY4/FLJ46536
Accession No.	Uniprot:Q6ZNE9Gene ID:285180
Uniprot	Q6ZNE9
GeneID	285180
Calculated MW	41
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only