

FLVCR2 Polyclonal Antibody FITC Conjugated

Catalog No: #C91183FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FLVCR2 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	FITC
Target Name	FLVCR2
Other Names	FLVCR2
Accession No.	Uniprot:Q9UPI3Gene ID:55640
Uniprot	Q9UPI3
GeneID	55640
Calculated MW	57
Concentration	1mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only