

## APPL1 Polyclonal Antibody FITC Conjugated

Catalog No: #C91192FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	APPL1 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/709
Conjugates	FITC
Target Name	APPL1
Other Names	APPL1
Accession No.	Uniprot:Q9UKG1 Gene ID:26060
Uniprot	Q9UKG1
GeneID	26060
Calculated MW	80
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only