

GTF2A2 Polyclonal Antibody FITC Conjugated

Catalog No: #C91193FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GTF2A2 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	51-109/109
Conjugates	FITC
Target Name	GTF2A2
Other Names	GTF2A2
Accession No.	Uniprot:P52657Gene ID:2958
Uniprot	P52657
GeneID	2958
Calculated MW	12
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only