

## B3GAT1 Polyclonal Antibody PE Conjugated

Catalog No: #C91195PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	B3GAT1 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	21-120/334
Conjugates	PE
Target Name	B3GAT1
Other Names	B3GAT1
Accession No.	Uniprot:Q9P2W7Gene ID:27087
Uniprot	Q9P2W7
GeneID	27087
Calculated MW	38
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only