

ZNF683 Polyclonal Antibody FITC Conjugated

Catalog No: #C91204FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZNF683 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu
Conjugates	FITC
Target Name	ZNF683
Other Names	ZNF683
Accession No.	Uniprot:Q8IZ20Gene ID:257101
Uniprot	Q8IZ20
GeneID	257101
Calculated MW	57
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only