

TNFR2 Polyclonal Antibody FITC Conjugated

Catalog No: #C91208FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TNFR2 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms
Immunogen Type	201-300/461
Conjugates	FITC
Target Name	TNFR2
Other Names	TNFR2
Accession No.	Uniprot:P20333Gene ID:7133
Uniprot	P20333
GeneID	7133
Calculated MW	46
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only