

INDOL1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91209Bio



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	INDOL1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-P IHC-F ICC IF
Species Reactivity	Hu
Immunogen Type	251-350/407
Conjugates	Bio
Target Name	INDOL1
Other Names	INDOL1
Accession No.	Uniprot:Q6ZQW0Gene ID:612129
Uniprot	Q6ZQW0
GeneID	612129
Calculated MW	45
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only