

Influenza A Polymerase basic protein 2 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91211Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Influenza A Polymerase basic protein 2 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC IF
Species Reactivity	Influenza A
Immunogen Type	681-759/759
Conjugates	Bio
Target Name	Influenza A Polymerase basic protein 2
Other Names	Influenza A Polymerase basic protein 2
Accession No.	Uniprot:PO3428Gene ID:956536
Uniprot	PO3428
GeneID	956536
Calculated MW	86
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only