

phospho-GRK2 (Ser670) Polyclonal Antibody FITC Conjugated



Catalog No: #C91214FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-GRK2 (Ser670) Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	PR(p-S)PV
Conjugates	FITC
Target Name	phospho-GRK2 (Ser670)
Other Names	phospho-GRK2 (Ser670)
Accession No.	Uniprot:P25098Gene ID:156
Uniprot	P25098
GeneID	156
Calculated MW	76
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only