

TRF4-2 Polyclonal Antibody FITC Conjugated

Catalog No: #C91218FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TRF4-2 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu
Immunogen Type	51-150/572
Conjugates	FITC
Target Name	TRF4-2
Other Names	TRF4-2
Accession No.	Uniprot:Q8NDF8Gene ID:64282
Uniprot	Q8NDF8
GeneID	64282
Calculated MW	63
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only