

KMT2C Polyclonal Antibody FITC Conjugated

Catalog No: #C91221FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	KMT2C Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	4701-4800/4911
Conjugates	FITC
Target Name	KMT2C
Other Names	KMT2C
Accession No.	Uniprot:Q8NEZ4Gene ID:58508
Uniprot	Q8NEZ4
GeneID	58508
Calculated MW	541
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200

Note: This product is for in vitro research use only