

KMO Polyclonal Antibody PE Conjugated

Catalog No: #C91223PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	KMO Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	51-150/486
Conjugates	PE
Target Name	KMO
Other Names	KMO
Accession No.	Uniprot:O15229Gene ID:8564
Uniprot	O15229
GeneID	8564
Calculated MW	56
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only