

KIAA0319L Polyclonal Antibody FITC Conjugated

Catalog No: #C91230FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	KIAA0319L Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	501-600/1049
Conjugates	FITC
Target Name	KIAA0319L
Other Names	KIAA0319L
Accession No.	Uniprot:Q8IZA0Gene ID:79932
Uniprot	Q8IZA0
GeneID	79932
Calculated MW	116
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only