

HPV16 E6 / HPV18 E6 Polyclonal Antibody Biotin Conjugated



Catalog No: #C91245Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HPV16 E6 / HPV18 E6 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-P IHC-F IF
Species Reactivity	Hu
Conjugates	Bio
Target Name	HPV16 E6 / HPV18 E6
Other Names	HPV16 E6 / HPV18 E6
Accession No.	Uniprot:P03126, P06463Gene ID:1489078 1489088
Uniprot	P03126, P06463
GeneID	1489078 1489088
Calculated MW	11
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only