

KAAG1 Polyclonal Antibody FITC Conjugated

Catalog No: #C91248FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	KAAG1 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu
Immunogen Type	31-84/84
Conjugates	FITC
Target Name	KAAG1
Other Names	KAAG1
Accession No.	Uniprot:Q9UBP8Gene ID:353219
Uniprot	Q9UBP8
GeneID	353219
Calculated MW	9
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only