

MASP1/MASP3 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91249Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MASP1/MASP3 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-F IF
Species Reactivity	Hu Ms Rat
Immunogen Type	20-120/699
Conjugates	Bio
Target Name	MASP1/MASP3
Other Names	MASP1/MASP3
Accession No.	Uniprot:P48740Gene ID:5648
Uniprot	P48740
GeneID	5648
Calculated MW	49/77
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only