

CX3CR1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91251Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CX3CR1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rat
Immunogen Type	151-250/355
Conjugates	Bio
Target Name	CX3CR1
Other Names	CX3CR1
Accession No.	Uniprot:P49238Gene ID:1524
Uniprot	P49238
GeneID	1524
Calculated MW	39
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200

Note: This product is for in vitro research use only