

phospho-MEK1 (Ser298) Polyclonal Antibody Biotin Conjugated



Catalog No: #C91255Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-MEK1 (Ser298) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F IF
Species Reactivity	Hu Ms Rat
Immunogen Type	PL(p-S)SY
Conjugates	Bio
Target Name	phospho-MEK1 (Ser298)
Other Names	phospho-MEK1 (Ser298)
Accession No.	Uniprot:Q02750Gene ID:5604
Uniprot	Q02750
GeneID	5604
Calculated MW	43
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only