

Melanoma gp100 Polyclonal Antibody Cy3 Conjugated

Catalog No: #C91263Cy3

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Melanoma gp100 Polyclonal Antibody Cy3 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu
Immunogen Type	21-120/661
Conjugates	Cy3
Target Name	Melanoma gp100
Other Names	Melanoma gp100
Accession No.	Uniprot:P40967Gene ID:6490
Uniprot	P40967
GeneID	6490
Calculated MW	68
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only