

DC-SIGNR1/CD209b Polyclonal Antibody FITC Conjugated

Catalog No: #C91265FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DC-SIGNR1/CD209b Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Ms
Immunogen Type	51-150/325
Conjugates	FITC
Target Name	DC-SIGNR1/CD209b
Other Names	DC-SIGNR1/CD209b
Accession No.	Uniprot:Q8CJ91Gene ID:69165
Uniprot	Q8CJ91
GeneID	69165
Calculated MW	37
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only