

Simian Rotavirus VP4 Polyclonal Antibody FITC Conjugated



Catalog No: #C91266FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Simian Rotavirus VP4 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Simian Rotavirus
Immunogen Type	201-300/776
Conjugates	FITC
Target Name	Simian Rotavirus VP4
Other Names	Simian Rotavirus VP4
Accession No.	Uniprot:P12473Gene ID:7011406
Uniprot	P12473
GeneID	7011406
Calculated MW	58/85
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only