

## EAAT2 Polyclonal Antibody PE Conjugated

Catalog No: #C91268PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	EAAT2 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/574
Conjugates	PE
Target Name	EAAT2
Other Names	EAAT2
Accession No.	Uniprot:P43004Gene ID:6506
Uniprot	P43004
GeneID	6506
Calculated MW	62
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only